

PRODUCT SPECIFICATION

SAFETYWARE 4-Ply Activated Carbon Face Mask (Earloop)

- This face mask provide nose and mouth protection against dust, pollen, bacteria and more.
- Made with premium quality Meltblown filter which offers high BFE & PFE of $\geq 99\%$.
- Outer layer is made of hydrophobic non-woven Polypropylene which provide good fluid splashes resistance while maintaining good breathability.
- The inner soft non-woven Polypropylene layer is hydrophilic which absorb droplets from the wearer effectively and makes the masks comfortable to wear.
- Build with activated carbon for filtration of organic vapours at nuisance level.
- Adjustable nose strip and flat elastic earloop to give basic protection to the user.
- Made by ultrasonic welding technology which provides strong bonding of mask & earloop.
- Comes with a smart pulls dispenser box which makes it convenient to pull up a mask.
- Features:
 - Soft & comfortable
 - Hypoallergenic
 - Low breathing resistance
 - Adaptable nose bar
 - Perfect fitting
 - UV sanitized
 - Disposable
- Made in Malaysia.

STORAGE

- The face mask shall maintain it's properties in a closed pack and stored in a cool dry condition.

Item Code: 1412EL-BLK



- SPLASH RESISTANT OUTER LAYER**
Hydrophobic - Protect Against Blood & Bodily Fluids
- ACTIVATED CARBON LAYER**
Filter Nuisance Odour
- MELTBLOWN FILTER LAYER**
Protect Against Bacteria & Particle
- SOFT INNER LAYER**
Hydrophilic - Absorb Droplets & Reduce Lens Fogging



PRODUCT SPECIFICATION

TECHNICAL DATA

COLOR	Black
STYLE	4-Ply with earloop
SIZE	17.5 cm x 9.5 cm
FILTRATION EFFICIENCY	BFE: $\geq 99\%$ PFE: $\geq 99\%$
TESTING METHOD	ASTM F2101: Bacterial Filtration Efficiency (BFE) ASTM F2299: Particle Filtration Efficiency (PFE)
SHELF LIFE	2 years from manufacture date
PACKAGING	50 pcs/ box, 30 boxes/ carton
CARTON DIMENSION	565 mm (L) x 415 mm (W) x 348 mm (H)

SAFETYWARE SDN. BHD. (632971-X)

Plot 237, Lengkok Perindustrian Bukit Minyak 3,
Bukit Minyak Industrial Estate,
14100 Simpang Ampat, Penang, Malaysia

Tel : +604-5023 882 | Fax : +604-5081 882
www.safetyware.com | info@safetyware.com

It is the product users' responsibility to conduct all the appropriate evaluations & testing to determine the suitability of SAFETYWARE products for use in a particular environment or purpose. SAFETYWARE disclaims all warranties other than as expressly provided.

Date created: 11/09/2020

Revision: 05

Date revised: 24/03/2025